


RAZ shower chair toilet compatibility

Standard WC- Adjusting the height of your RAZ for clearance

The table below shows the seat height setting you will need to clear the listed toilet bowl heights. It is important to make a note of the toilet bowl height when carrying out a client assessment and submitting this information when getting a quotation or ordering from core stock.


| Toilet bowl height (mm) | 380 | 400 | 420 | 440 | 460 | 480 | 500 | 520 |
|-----------------------------|-----|-----|-----|-----|-----|-----|-----|-----|
| AT standard seat height | 530 | 530 | 560 | 560 | 580 | 610 | 640 | 640 |
| ATP standard seat height | 530 | 530 | 560 | 560 | 580 | 610 | 640 | 640 |
| AP/SP standard seat height | 500 | 520 | 540 | 560 | 580 | 600 | | |
| AT600 5" castor seat height | 570 | 570 | 570 | 570 | 570 | 600 | | |
| AP600 5" castor seat height | 490 | 510 | 540 | 570 | 570 | 600 | | |
| SP600 5" castor seat height | 490 | 510 | 540 | 570 | 570 | 600 | | |



Measurement from floor to bottom of seat

Lining up your RAZ with the toilet bowl

Getting this measurement correct is important so that the seat aperture lines up with the toilet opening. Remember to consider waste pipe location and the toilets distance from the wall.



| Rearmost chair point to rear of seat aperture | |
|---|-----|
| AT (attendant tilt) | 290 |
| ATW (attendant tilt wide) | 290 |
| ATP (attendant tilt paediatric) | 230 |
| AP (attendant propel) | 180 |
| AP Auto (attendant propel wide) | 180 |
| SP (self propel) | 280 |
| AT600 (attendant tilt bariatric) | 350 |
| AP600 (attendant propel bariatric) | 255 |
| SP600 (self propel bariatric) | 380 |

Best Practice

A common requirement is to have the skirt accessory added to the underside of the seat and the option of removing the pan sliders to provide a closer fit with the toilet. The standard moulded seat shaping contributes to a much improved wash experience.



Bio Bidet

The RAZ range is not compatible with a Bio Bidet.



RAZ shower chair toilet compatibility

Auto WC Compatibility

See below which RAZ chairs are compatible with popular automatic toilets. If your clients automatic toilet is not listed below then please contact our sales team to discuss.

AT



✓ Geberit Mera

ATW



- ✓ Geberit Mera
- ✓ Geberit 7000
- ✓ Aquaclean 8000
- ✓ Closomat Samoa / Palma Vita

AP & SP



✓ Geberit Mera

AP GEB



- ✓ Geberit Mera
- ✓ Geberit 7000
- ✓ Aquaclean 8000
- ✓ Closomat Samoa / Palma Vita

AT/AP/SP 600



✓ Geberit Mera

RAZ ATP



✓ Geberit Mera



Assessment Considerations

- The size and shape of the room
- Turning circle
- Obstructions around the toilet- grab handles, basins etc
- Soil stack location and direction
- Flooring and thresholds
- Remember the RAZ will need to be in an upright position when over a WC.

Contact Wealden Rehab for assessment forms